


Oral/Max


Dedicated Saws



Expansive Drilling Capabilities



MicroDrill and Attachments



U2 Drill with Elite Attachments

Integrated Irrigation



- Specialty-driven handpieces
- Increased power, precision cutting
- Complete OR standardization

Impaction Drill

- 85,000 RPM
- 83% more torque
- Electronic Torque Feedback for optimal speed and performance




Universal Driver

- Trigger Oscillating Feature
- 40% More Torque
- 200% More Power



Oscillating Saw

- 32% Lighter, 15% Shorter for sleek, ergonomic design
- Tool-less crescentric, intra-oral blades
- New blade mount



Sagittal Saw

- Rotating, Indexing Head
- 40% More Power than TPS Sagittal Saw
- New Mount to swing large blades (Rayhack)



Reciprocating Saw

- 23% Lighter, 11% Shorter for a more ergonomic design



CORE Oral/Max Surgery Ordering Information

PART #	DESCRIPTION	PART #	DESCRIPTION
5400-300	CORE Impaction Drill	5400-099-000	CORE Universal Driver
2296-301-000	Impaction Bur Guard	4100-120-000	T-Latch Attachment
2296-300-000	Impaction Bur Shield	4100-125-000	Pin Collet (.078" - .125" Diameter)
5400-15	CORE Micro Drill	4100-062-000	Wire Collet (.028"-.071" Diameter)
5400-9	CORE Hand Switch	4100-110-000	AO Small Drill
5100-15-43	MD Series 4:1 Straight Reducer	4100-131-000	1/4" Drill Chuck (w/ Key)
5100-15-250	MD Series Medium Straight Attachment	4100-132-000	5/32" Drill Chuck
5100-15-252	MD Series Medium Angled Attachment	4100-133-000	1/4" Keyless Drill Chuck
5100-15-270	MD Series Long Straight Attachment	4100-135-000	Hudson®/Modified Trinkle Drill
5100-15-272	MD Series Long Angled Attachment	4100-160-000	Trinkle Drill
5400-34	CORE Sagittal Saw	4100-210-000	AO Large Reamer
5400-31	CORE Oscillating Saw	4100-225-000	Pin Collet (.126"-.157" Diameter)
5400-37	CORE Reciprocating Saw	4100-231-000	1/4" Reamer Chuck (w/ Key)
		4100-235-000	Hudson®/Modified Trinkle Reamer
		4100-260-000	Trinkle Reamer
		4100-260-000s4	DHS/DCS quick lock
		4100-400-000	Sagittal Saw Attachment
		4100-700-000	High Speed Bur Attachment



For complete CORE System ordering information, contact your local Stryker Sales Representative.

stryker®
MedSurgGroup

See Bur and Blade literature for complete listing of available cutting accessories

Instruments Division
4100 East Milham Ave
Kalamazoo, MI 49001
(800) 253-3210

Endoscopy Division
5900 Optical Court
San Jose, CA 95138
(800) 624-4422

Leibinger Division
4100 East Milham Ave
Kalamazoo, MI 49001
(800) 962-6558